




Sacral Neuromodulation with the InterStim[®] System


PATIENT MANAGEMENT GUIDE



Innovating for life.

This document offers suggested best practices for implementing Manage-By-Fact in a practice. Programming decisions are the responsibility of the clinician and are based on the ability to establish settings that will provide optimal patient symptom relief, minimize patient discomfort and maintain neurostimulator battery life to the best possible extent.





Patients with a Patient Management Worksheet on file:

- Instruct patient to record symptoms using a Voiding Diary for a minimum of 3 days and call back with results.
 - If >50% improvement from baseline: Review patient expectations. Reassure patient that stimulation is on. Turning device off may indicate therapy is working as intended.
 - If <50% improvement from baseline: Instruct patient to try the other program(s) in iCon™ Programmer.
- If necessary, repeat step 1 to evaluate other programs in the iCon™ Programmer. Advise patient to evaluate each program for a minimum of 3 days and record symptoms using a Voiding Diary.
- If all 4 programs in iCon™ Programmer have been tried, schedule an office visit.

At appointment:

- Review Voiding Diary information and fill out **Patient Management Worksheet**.
- Run impedance check & review Therapy Usage using the N'Vision® Clinician Programmer.
- Load other programs using the N'Vision® Clinician Programmer.
- Advise patient to activate 1st program in iCon™ Programmer for a minimum of 3 days.

Post-appointment, instruct patient to:

- If necessary, repeat with trying the other program(s) in iCon™ Programmer for a minimum of 3 days and record symptoms using a Voiding Diary.
- Schedule an office visit if necessary.

Physicians can make the decision to REPROGRAM, REVISE or REMOVE the neurostimulator at any time.

Patients without a Patient Management Worksheet on file:

Instruct patient to:

- Turn device off with iCon™ Patient Programmer (if device is on).
- Record symptoms using a Voiding Diary for a minimum of 3 days.
- Schedule an office visit.

At appointment:

- Collect baseline diary information and fill out **Patient Management Worksheet**.
- Run impedance check & review Therapy Usage using the N'Vision Clinician® Programmer.
- Load 4 programs into iCon™ Programmer.
- Advise patient to activate 1st program in iCon™ Programmer for a minimum of 3 days.

Post-appointment, instruct patient to:

- If necessary, repeat with trying the other program(s) in iCon™ Programmer for a minimum of 3 days and record symptoms using a Voiding Diary.
- Schedule an office visit if necessary.

Physicians can make the decision to REPROGRAM, REVISE or REMOVE the neurostimulator at any time.

Indications for Use:

InterStim® Therapy for Urinary Control is indicated for the treatment of urinary retention and the symptoms of overactive bladder, including urinary urge incontinence and significant symptoms of urgency-frequency alone or in combination, in patients who have failed or could not tolerate more conservative treatments.

The following Warning applies only to InterStim Therapy for Urinary Control:

Warning: This therapy is not intended for patients with mechanical obstruction such as benign prostatic hypertrophy, cancer, or urethral stricture.

InterStim® Therapy for Bowel Control is indicated for the treatment of chronic fecal incontinence in patients who have failed or are not candidates for more conservative treatments.

Contraindications for Urinary Control and for Bowel Control: Diathermy. Patients who have not demonstrated an appropriate response to test stimulation or are unable to operate the neurostimulator.

Precautions/Adverse Events:

For Urinary Control: Safety and effectiveness have not been established for bilateral stimulation; pregnancy, unborn fetus, and delivery; pediatric use under the age of 16; or for patients with neurological disease origins such as multiple sclerosis.

For Bowel Control: Safety and effectiveness have not been established for bilateral stimulation; pregnancy, unborn fetus, and delivery; pediatric use under the age of 18; or for patients with progressive, systemic neurological diseases.

For Urinary Control and for Bowel Control: The system may be affected by or adversely affect cardiac devices, electrocautery, defibrillators, ultrasonic equipment, radiation therapy, MRI, theft detectors/ screening devices. Adverse events include pain at the implant sites, new pain, lead migration, infection, technical or device problems, adverse change in bowel or voiding function, and undesirable stimulation or sensations, including jolting or shock sensations. For full prescribing information, please call Medtronic at 1-800-328-0810 and/or consult Medtronic's website at www.medtronic.com. Product technical manual must be reviewed prior to use for detailed disclosure.

USA Rx Only. Rev 0409

www.medtronic.com

United States of America

Medtronic Neuromodulation
710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA
Tel. +1-763-505-5000

Asia-Pacific

Medtronic International, Ltd.
Suite 1106-11, 11/F, Tower 1, The Gateway
25 Canton Road, Tsimshatsui
Kowloon
Hong Kong
Tel. +852-2919-1362

Canada

Medtronic of Canada Ltd.
99 Hereford Street
Brampton
Ontario L6Y 0R3
Canada
Tel. +1-905-460-3800

Europe

Medtronic International Trading Sàrl
Route du Molliou 31
Case Postale 84
CH-1131 Tolochenaz
Switzerland
Tel. +41-21-802-7000

Australia

Medtronic Australasia Pty. Ltd.
97 Waterloo Road
North Ryde, NSW 2113
Australia
Tel. +61-2-9857-9000

