



LIGHTSHEER

This is how confidence feels

LightSheer[®] Laser Hair Removal

A collection of Before & After images
from around the world



We don't just build technology, we build confidence.

What we do changes the way people feel about themselves.
We don't just remove unwanted hair, we remove doubts and worries
and the barrier between people and what they can be.
We create, design and manufacture confidence.

This is how **confidence** feels



Laser Hair Removal 805nm

Lower Abdomen



Baseline



Followup after 6 treatments

LightSheer DUET, High Speed handpiece, 805 nm wavelength Skin type III, Fluence: 8.5-10J/cm², Pulse duration: 30 ms

Courtesy of Daniela Cangelosi, MD

Back and Shoulders



Baseline



After 2 treatments

LightSheer DUET, High Speed handpiece Fluence: 10J/cm², Pulse duration range: 30-70 ms (Custom mode), 805 nm wavelength

Courtesy of Mitchel Goldman, M.D.

Axilla



Baseline



6 months follow up after 3 treatments

LightSheer DUET, High-Speed handpiece Skin type IV, Brown hair , Fluence: 11 J/cm², Pulse duration: 60-66ms, 805 nm wavelength

Courtesy of Suzanne Kilmer, M.D.

Laser Hair Removal 805nm

Lower Abdomen



Baseline



6 months follow up after 5 treatments

LightSheer DUET, High Speed handpiece, 805 nm wavelength

Courtesy of Girish Munavalli, M.D., M.H.S

Chest



Baseline



2 months follow-up after 5 treatments

LightSheer DESIRE, High Speed handpiece, 805nm wavelength, Fluence: 7.5-8 J/cm², Pulse Duration: 100 ms

Courtesy of Shweta Virmani, M.D.

Upper Back and Shoulders



Baseline



6 months follow up after 5 treatments

LightSheer INFINITY, High Speed handpiece, 1060 nm wavelength Skin type IV, Fluence: 6-7.5 J/cm², Pulse Duration: 30 ms

Courtesy of Sushil Tahiliani, M.D.

Axilla



Baseline

LightSheer INFINITY, High Speed handpiece, 1060 nm wavelength



6 months follow up after 5 treatments

LightSheer INFINITY, High Speed handpiece, 1060 nm wavelength Skin type IV, Fluence: 12 J/cm², Pulse Duration: 50-100 ms

Courtesy of Sushil Tahiliani, M.D.

Axilla



Baseline



6 months follow up after 5 treatments

LightSheer INFINITY, High Speed handpiece, 1060 nm wavelength Skin type IV, Fluence: 12 J/cm², Pulse Duration: 50-100 ms

Courtesy of Sushil Tahiliani, M.D.

Lower Arms



Baseline



6 months follow up after 5 treatments

LightSheer INFINITY, High Speed handpiece, 1060nm wavelength Skin type V, Fluence: 7.5-12.5 J/cm², Pulse Duration: 30-50 ms

Courtesy of Sushil Tahiliani, M.D.