



Hand to Shoulder Surgery & Reconstructive Microsurgery • [www.achssurgeons.com](http://www.achssurgeons.com)  
Phone (602) 258-4788 • Fax (602) 258-5131

# Tips for Optimal Placement of Immobilizer

## Vivek Agrawal, M.D.

### Harness Position at Opposite Shoulder

---



## Abductor Pillow Placed High in Axilla

---

- Valve Access Panel should be positioned posteriorly with BREG logo facing forward. Adjust harness straps to suspend pillow and allow operative arm to rest on pillow tension free.
- Arm and wrist Velcro straps are then positioned and applied loosely for support.



## Pillow Positioned to Keep Elbow Forward of Midline

---

- Once the pillow is placed high in axillar the arm should be draped loosely over pillow for comfort.
- Elbow should remain forward of midline when adjusting and applying the arm and wrist holders.



## Posterior View

---

- Larger Velcro loop is used for arm to help prevent it from falling backwards beyond the midline.



## Side View

---

- Ball may be placed to allow wrist and hand to rest comfortably.

