

ROSENBAUM POCKET VISION SCREENER

95

distance
equivalent
 $\frac{20}{800}$

874

Point
Jaeger $\frac{20}{400}$

2843

26 16 $\frac{20}{200}$

638 E W E X O O

14 10 $\frac{20}{100}$

8 7 4 5 E M W O X O

10 7 $\frac{20}{70}$

6 3 9 2 5 M E E X O X

8 5 $\frac{20}{50}$

4 2 8 3 6 5 W E M O X O

6 3 $\frac{20}{40}$

3 7 4 2 5 8 E W E X X O

5 2 $\frac{20}{30}$

9 3 7 8 2 6 W M E X O O

4 1 $\frac{20}{25}$

4 2 8 7 3 9 E M W O O X

3 1+ $\frac{20}{20}$

Card is held in good light 14 inches from eye. Record vision for each eye separately with and without glasses. Presbyopic patients should read thru bifocal segment. Check myopes with glasses only.

DESIGN COURTESY J.G. ROSENBAUM, M.D.

PUPIL GAUGE (mm.)

