

# What to Expect with Theradome Laser Hair Therapy

- 1 Slow down hair loss (4-12 weeks)
- 2 Thicken existing hair (12-24 weeks)
- 3 Renew and grow new hair (24 - 36 weeks)



*Less Hair  
In Your  
Brush*



*Hello  
Hairless  
Pillow*



*No More  
Hair In The  
Drain*

4 out of 5 hair loss professionals recommend  
**AWARD WINNING**  
Theradome



Find out how Theradome can restore your hair and your confidence

Keep the Hair You Have and **Grow New Hair**



With the world's most **technologically advanced** laser phototherapy device

Learn more at: [www.theradome.com](http://www.theradome.com)

**theradome**  
HAIR GROWTH TECHNOLOGY

# Why Choose Theradome?



Theradome is the leading laser phototherapy (LPT) device for hair loss and is the go-to choice for patients and physicians.

- ✓ Fast results vs. other LPT devices
- ✓ FDA cleared for both women and men
- ✓ More Androgenetic Alopecia scalp coverage than any other device
- ✓ Convenient at home use
- ✓ Easy one touch activation
- ✓ Hands-free, cordless and lightweight
- ✓ Voice guided treatment counter
- ✓ One charge up to 7 treatments
- ✓ Heat management that enables maximum dose for effective treatment, safely and in comfort



## Life-Changing Results

What real customers have achieved with Theradome laser hair restoration devices\*



DEBRA - 52 years old



FERNANDO - 48 years old



CELESTE - 48 years old



NANCY - 65 years old



WARREN - 71 years old

\*Individual results may vary

## How Does Laser Phototherapy Work?

### Theradome's targeted lasers


Deliver a precise dose of energy to your hair follicles, bio-stimulating your mitochondria (the "power plants" of your cells) in order to stop hair loss, thicken hair follicles, and grow new hair.

### A precise wavelength proven to get results


Theradome helmets are specifically engineered to deliver laser light at 680 nm wavelengths, the gold standard in stimulating hair growth, according to the Swedish Laser Medical Society.

Most other LPT devices do not meet this same standard and some don't even use real lasers at all, relying on cheap LEDs instead.

### What's the difference between LASERS and LEDs?



**Medical Grade LASERS**  
Theradome medical grade LASERS penetrate the scalp up to 5mm for maximum follicle stimulation, more than any other at home laser hair growth device.



**LEDs**  
Some laser hair loss products use LEDs. LED beams diverge without reaching the base of hair follicles.