

The results obtained with FibroScan, along with other diagnostic tests, will help your physician monitor your liver health.



To learn more  
Ask your physician  
or hepatologist.

## My FibroScan examination

### APPOINTMENT

Date: \_\_/\_\_/----

Time: \_\_/\_\_/----

Place: \_\_\_\_\_



Fast,  
painless  
assessment  
of your  
liver health



Sandhill Scientific, Inc  
9150 Commerce Center Circle, #500  
Highlands Ranch, CO 80129 USA  
800.468.4556 303.470.7020  
Fax: 303.470.2975

[FibroScan502touch.com](http://FibroScan502touch.com)

**FibroScan** touch **502**




*«Fast & painless, FibroScan makes the monitoring of my disease more convenient.»*

## What is FibroScan?

An examination with FibroScan, which uses a technology called Vibration Controlled Transient Elastography (VCTE™), is a non-invasive test used to assess liver stiffness.

FibroScan is indicated for noninvasive measurement of shear wave speed and equivalent stiffness in the liver. The shear wave speed or stiffness may be used as an aid to clinical management of your liver disease.

The result is immediate and shows the evolution of the health of your liver.



FIBROSCAN MEASURES LIVER STIFFNESS.

## Is it painful? No.

The FibroScan examination is painless, quick and easy. During measurement, you feel a slight vibration on the skin at the tip of the probe.

## What does the FibroScan examination consist of?

- You lie on your back, with your right arm raised behind your head. Your technician or nurse applies a water-based gel to the skin and places the probe against your skin with a slight pressure.
- The examination includes 10 consecutive measurements made at the same location.
- The result is delivered at the end of the examination. Your doctor will interpret the result.

## Can I eat prior to the test?

You should not eat within 2 or 3 hours of the FibroScan test.

## What do I wear?

The FibroScan test is performed on your right side under your arm at the level of your breast bone. Patients should wear a blouse or shirt which can be partially unbuttoned to expose their right side above the waist and below the breast bone.

## How long is the test?

The entire FibroScan test session will require 15-20 minutes and provides immediate results.

*« My physician suggested that I try the FibroScan »*



*FibroScan probes are adapted to all body-types.*



## What to do with the result?

Your physician interprets the result according to your history and your disease.

## What difference does FibroScan make for me?

- FibroScan examination is easy and fast.
- The exam is painless and non-invasive.
- The FibroScan test can be safely repeated.