

Great-looking skin doesn't have to be complicated

With all the invasive procedures, injectables, chemicals, and countless lotions and potions out there, isn't it time for a simple solution that works?

The TempSure Envi treatment minimizes facial fine lines and wrinkles, tightens* skin and improves the appearance of cellulite. It's safe for all skin types and levels of sun exposure, so it fits perfectly into any beauty maintenance routine, any time of year.

With TempSure Envi, there's no surgery, no needles, and no downtime, so you can immediately get back to what you were doing, and look great doing it.

Welcome to a simple way to maintain beautiful skin. Welcome to TempSure Envi.

Ask your provider about
TempSure Envi today!
For more information, visit
www.beautifullysimple.com

BEAUTIFUL SKIN

BEAUTIFULLY Simple



TempSure *Envi*

The TempSure Generator has the following indications for use:
The 10mm, 15mm, and 20mm TempSure Envi handpieces are indicated for nonablative treatment of mild to moderate facial wrinkles and rhytids. The 18mm, 25mm, and 30mm TempSure Envi handpieces provide heating for the purpose of elevating tissue temperature for selected medical conditions such as temporary relief of pain, muscle spasms, and increase in local circulation.
The TempSure Envi massage device is intended to provide a temporary reduction in the appearance of cellulite.

CYNOSURE
A Hologic Company

*Through soft tissue coagulation ©2017 Hologic, Inc. All Rights Reserved. Cynosure is a registered trademark of Cynosure, Inc. TempSure is a trademark of Cynosure, Inc. Hologic owns exclusive rights to the proprietary use of sonography without written permission of Hologic is prohibited. Model for skin shown is unretouched only. Results may vary. 021-7027-020 Rev 3/1918



How does it work?

The TempSure Envi treatment delivers radiofrequency energy to gently heat your skin. Increasing the skin's temperature for a defined, therapeutic time will trigger a natural response for your skin to create new collagen. The new collagen fibers are tight*, dense, and neatly organized giving you tighter*, younger-looking skin.

What are your problem areas?

Application areas include forehead lines, frown lines, crow's feet, smile lines, wrinkles in the chin, and cellulite on the body.



FAQs

Does it hurt?

The TempSure Envi treatment has been carefully designed with comfort in mind. 99% of patients described the treatment as comfortable*.

How long are treatments?

Treatments last about 20-60 minutes depending on the area being treated. There is absolutely no downtime allowing you to return to your normal activities immediately.

When will I see results?

Instantly! You will leave with a natural glow, and what patients describe as a "tighter" feeling*. After the treatment your collagen will continue to rebuild, enhancing your skin's appearance overtime.

Is this right for me?

Treatments are extremely safe and can be performed on all skin types. Unlike other technologies, treatments can be performed all year round.

Why TempSure Envi?

The TempSure Envi treatment is a simple solution for beautiful skin. It provides a soothing patient experience, while delivering beautiful results. Each treatment delivers consistent heat, giving you an effective, yet comfortable treatment every time!

How many treatments are needed?

A series of treatments is always recommended, but every patient is different. The optimum number of treatments can vary, so please work with your provider to develop a plan that is right for you!

*Through soft tissue coagulation
In-house study

NO SURGERY
NO INJECTIONS
NO DOWNTIME
REAL RESULTS



Before, and one month after 2nd treatment. Courtesy of S. Doherty, MD.